**sartorius stedim**  
biotech

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tech.pubs@sartorius-stedim.com  
www.sartorius-stedim.com

Systemvoraussetzungen:

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The enclosed CD contains the operating instructions for SARTOFLOW® Alpha plus SU as a PDF file in German, English, French, Spanish and Italian.

Further languages on request.

If the CD is missing, you can obtain a copy form us by specifying the order number:

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If you require further information on crossflow filtration, other than mentioned in this User Manual, please contact us.

Reproduction, translation and duplication in any form, in part or in full, requires the written permission of the Sartorius Stedim Biotech GmbH.

For technical information write to:

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
D-37079 Goettingen  
Internet: [www.sartorius-stedim.com](http://www.sartorius-stedim.com)

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Version: Juni 2010  
Bestellnummer: 85032-537-29

# Contents

<b>1. Structure and Application</b>	<b>10</b>	<b>3. Operation of the Filtration System</b>	<b>23</b>
1.1 General Notes	10	3.1 Starting the Unit   Emergency Stop Function	23
1.2 System Structure	10	3.2 Display	23
<b>2. Installation and Starting Up</b>	<b>12</b>	3.2.1 Touch Fields of the Start Screen	23
2.1 In General	12	3.3 Modes of Operation	24
2.1.1 Assembly Surface	12	3.3.1 Manual Operation	24
2.1.2 Preparation of the Workplace	13	3.3.2 Semi-Automatic Operation	24
2.2 General Connections and Installation Material of the Filtration Module	13	3.3.3 Menu Bar	24
2.2.1 Mains Supply	13	<b>4. Appendix</b>	<b>41</b>
2.2.2 Communication (Ethernet)	13	4.1 Technical Data	41
2.2.3 Sewage Discharge	13	4.2 Cleaning and Care	41
2.2.4 Discharged Air   Ventilation	13	4.3 Maintenance and Service	41
2.3 DCU4 Tower	14	4.4 Causes of Malfunction	42
2.3.1 Structure of the DCU4 Tower	14	4.5 Disassembly and Disposal	43
2.3.2 Installation and Alignment of DCU4-Tower	14	4.6 Warranty	44
2.3.3 Signal and Power Supply of the DCU4 Tower	15	4.7 EC Declaration of Conformity	44
2.4 Filtration Unit	15	4.8 Declaration of Decontamination	44
2.4.1 Structure of the Filtration System	15	4.9 Contact Addresses	44
2.4.2 Assembly and Alignment of Filtration Unit	16	4.10 Information and Instructions on Disposal	45
2.4.3 Inspection of the filtration unit	16		
2.4.4 Signal and Power Supply of the Filtration Unit	17		
2.4.5 Mounting of the Recirculation Container (Bag as Recirculation Container)	17		
2.4.6 Installation of the presterilized Single use Crossflow Loop	18		

# Introduction

This operating manual describes the installation and commissioning, the connection of the energy | supply modules as well as the operation of a SARTOFLOW® filtration system.

Information on the different components you can find in the external operating manuals.

## Application of the SARTOFLOW® Filtration System

The SARTOFLOW® filtration system is suited for the continuous filtration of suspensions and liquid media within the laboratory and production scale. By using different kinds of Self Contained Units the filtration system can be used for microfiltration as well as for ultrafiltration purposes.

For special processes the use may be permitted only under certain circumstances. As concerns e. g. the biological safety there are special requirements on the workplace, the instrumentation, the handling of the components, the safety of the personnel and the work environment. This operating manual does not specify such requirements and other regulations which are obligatory by law or otherwise.

The operating of the filtration system requires the usual knowledge of an expert for the handling of micro organisms, cells and contaminated items. If there are special risks for users and the working environment at the handling of components, the operating manual contains danger, warning, caution and safety notes. They apply in addition to the regulations for the use at the workplace.

The SARTOFLOW® filtration system has been designed for stationary operation purposes only.

## Signs and symbols

**Danger notes bear the preceding symbol and they are marked in red colour. If you disregard these danger notes damages of the equipment or other material damages or personal injury with a high risk may result.**

**Notes of warning bear the preceding symbol and they are marked in orange colour. If you disregard these notes of warning damages of the equipment or other material damages and personal injury with a medium risk may result.**

**Caution notes bear the preceding symbol and they are marked in yellow colour. If you disregard these caution notes damages of the equipment or other material damages and personal injury with a low risk may result.**

**Safety notes and general notes bear the preceding symbol. It identifies all steps which have to be carried out with extreme care. It appears if special aspects have to be considered.**

The refers to any information included in this or in other documents. The brackets specify the sources of figures, paragraphs or documents.

These paragraphs contain operating steps which have to be carried out one after another.

## Introduction

The SARTOFLOW® microfiltration or ultrafiltration system is a part of the product program of Sartorius Stedim.  
If you have any questions concerning this system or the rest of the product program of Sartorius Stedim you can contact the following address.

Sartorius Stedim Systems GmbH  
Robert-Bosch-Str. 5-7  
D - 34302 Guxhagen, Germany  
Phone: +49.5665.407.0  
Fax: +49.5665.407.2201

## Notes to this documentation

All information on the SARTOFLOW® filtration system refers to the time of the publication of this document. We reserve the right to modify different components of the filtration system compared with the information and modifications of this document without separate notice. Third parties are not allowed to alter the operating manual in some way without our previous written authorization. The same applies to the copying or the use of the document for other purposes.

# Safety notes

## General notes

Please note:

- The operators must be technically qualified for the field of application and they must be familiar with the handling of the unit and the resulting dangers.
- Your company must hang up safety regulations on existing dangers. They must be visible (at the application) for the employees and the work environment at the place of assembly. If necessary suitable safety devices have to be provided.
- It is generally necessary to wear suited work wear and personal protective equipment such as gloves, protective goggles and eventually a breathing protection.
- Unauthorized employees are not allowed to work at the unit.

## Notes for the save use of the filtration system

Please note:

- At the work at and with the SARTOFLOW® filtration system you only have to use such devices and peripheral devices which have been released for this kind of use by Sartorius Stedim.

### Risk!

**Be very careful when handling with the recirculation bag. You may injure or unintentionally spill the media.**

---

- If modifications are carried out at the filtration system you have to handle the recirculation system with care. Please check the proper condition of all parts before using them. Damaged parts must not be installed.  
Optionally mounted parts, external connections and components have to be fastened carefully and protected from loosening unintentionally.
- Please take care of the weight and the required space of the components.  
The installation surface must be able to carry the completely equipped filtration system ► „see chapter 2.1.1 | 2.1.2“. You have to use suitable transport aids for the transport to the installation location and for relocations.
- Ensure that the supply voltage as well as any other supply powers correspond with the device specifications ► „see chapter 2.3.3, chapter 4.1“.

Please note:

---

**Risk of breakage or bursting!**

**You must charge the recirculation bag only with the maximum pressure of 300 mbar.**

**At all plastic or similar materials (e. g. bags) it is not recommended to use pointy items (scissors, knife, screwdriver etc.) since they may break.**

---

- The SARTOFLOW® filtration system must only be operated with the admissible energies (power supply) and the provided performance data.
- The technical data are described in a separate chapter
  - ▶ „see chapter 4.1“ of the operating manual.
- You must always control the processes if critical situations may occur. You must make arrangements in order to avoid any dangers or damages due to process disturbances.
- ▶ „After finishing the filtration process and before any maintenance works you have to observe the safety regulations for contaminated devices“.

---

**Interventions and technical modifications at the DCU4 tower, the control modules, the supply modules, the drive modules, the vessels and the peripheral devices as well as at all adapters for the connection to the peripheral devices are not allowed if you have not been authorized by Sartorius Stedim Systems GmbH.**

---

Safety notes

- Only authorized persons are allowed to carry out service and repair works. You must only use spare parts which have been released for the existing SARTOFLOW® filtration system by Sartorius Stedim Systems GmbH.

## 1. Structure and Application

### 1.1 General Notes

The SARTOFLOW® Alpha plus SU filtration system has been designed for Single use applications as well as for the use at the production of small volumes. The automatically mode of operation of the SARTOFLOW® Alpha plus SU filtration system guarantees repeatable results.

SARTOFLOW® Alpha plus SU can be equipped with a gamma presterilized loop. All product-contact surfaces like self contained UF|MF unit, pressure domes, flow meters, valves, bags, sterile connectors and tubing support a maximum of flexibility and process security.

The SARTOFLOW® Alpha plus SU is particularly suitable for contract manufacturers, research and development applications, small scale production cycles and companies with frequently changing productions.

### 1.2 System Structure

The SARTOFLOW® Alpha plus SU consists of the following:

No.	Modules	Description
1	SARTOFLOW® Alpha plus SU filtration unit with standard filter holder	Filtration unit with standard instrumentation, and bag stand. Manual filter holder
2	SARTOFLOW® Alpha plus SU filtration unit with hydraulic filter holder	Filtration unit with standard instrumentation, and bag stand. Hydraulic filter holder
3	DCU4 Tower 230 VAC   60 Hz	
4	DCU4 Tower 120 VAC   50 Hz	

#### Consumables

Order number	Description	
SFA-SU-1463901	UF 10 kDa PESU, 0.1 m <sup>2</sup>	components of bag assembly: 10 l recirculation bag, pressure domes, flow meters, valve bodies, OPTA connectors and tubing
SFA-SU-1463921	UF 10 kDa PESU, 0.2 m <sup>2</sup>	
SFA-SU-1463931	UF 10 kDa PESU, 0.3 m <sup>2</sup>	
SFA-SU-1465901	UF 30 kDa PESU, 0.1 m <sup>2</sup>	
SFA-SU-1465921	UF 30 kDa PESU, 0.2 m <sup>2</sup>	
SFA-SU-1465931	UF 30 kDa PESU, 0.3 m <sup>2</sup>	
TBD	MF 0.1 µm PESU, 0.1 m <sup>2</sup>	
TBD	MF 0.1 µm PESU, 0.2 m <sup>2</sup>	
TBD	MF 0.1 µm PESU, 0.3 m <sup>2</sup>	

Figure 1-1: SARTOFLOW® Alpha plus SU

1. Hardware: DCU4 control tower, Filtration unit with standard instrumentation, filter holder and bag stand.
2. Gamma presterilized Loop consisting of Self Contained Unit MF or UF, recirculation bag, pressure domes, flow meters, valves, OPTA sterile connectors and tubing.

## 2. Installation and Starting Up

### 2.1 In General

## 2. Installation and Starting Up

The SARTOFLOW® Alpha plus SU filtration unit is supplied to the customer after a thorough functional check.

In case the SARTOFLOW® Alpha plus SU filtration unit should have been damaged during the transport or after the delivery and installation you have to pay attention to the following details due to eventual warranty claims (and also due to other defects).

- Check the completeness of the supplied material according to your order.
- Thoroughly check all components for damages.
- Report the missing parts or any transport damages to your local representation of Sartorius Stedim Systems GmbH or to Sartorius AG.
- Please check the operability of the SARTOFLOW® Alpha plus SU filtration unit before starting the first filtration process. Please thoroughly record eventual malfunctions. A record may facilitate you and Sartorius the handling of eventual claims.

#### 2.1.1 Assembly Surface

The lower figure shows the SARTOFLOW® Alpha plus SU as complete unit with the minimum external dimensions.

All dimensions in mm.

Figure 2-1: External dimensions of the SARTOFLOW® Alpha plus SU

### 2.1.2 Preparation of the Workplace

The installation location must be suited for the filtration system. The installation location must offer enough space and the bottom must be able to carry the weight of the completely equipped unit.

- When installing the SARTOFLOW® Alpha plus SU you must take care for sufficient space for the connection to the laboratory devices as well as other necessary modules.

Check the correct dimensioning of the following connections:

- Mains voltage.

## 2.2 General Connections and Installation Material of the Filtration Module

### 2.2.1 Mains Supply

#### **Voltage damage!**

**The connection to a wrong power supply may lead to damage of the unit. Before starting the device the correct power supply has to be checked carefully.**

**Check at the type plates of the systems components supplied to you if the voltage type is correct and if the power cables are provided with the right plug necessary for your kind of use.**

If you have got a wrong voltage type or a power cable with wrong plugs, please contact your customer service at Sartorius Stedim Systems GmbH.

**The DCU4 tower of the SARTOFLOW® Alpha plus SU is available in two voltage types.**

**230 volts (±10 %), 50 | 60 Hz**

**or**

**110 volts (±10 %), 50 | 60 Hz**

### 2.2.2 Communication (Ethernet)

In order to communicate with superior central computer systems (SCADA) the SARTOFLOW® Alpha plus SU has got an Ethernet connection in the DCU4 tower module.

### 2.2.3 Sewage Discharge

The discharge of sewage is carried out by means of respective connections (part of presterilized Loop).

### 2.2.4 Discharged Air | Ventilation

The pressure equalization within the filtration unit is carried out by means of a ventilation device (part of presterilized Loop).

## 2.3 DCU4 Tower

### 2.3.1 Structure of the DCU4 Tower

The „DCU4 tower“ serves as digital measurement and control system which is connected to the filtration unit.

► „Technical description: see external operating manual „DCU4 Tower“.


Figure 2-2: Front view of the DCU4 tower

1. Touch panel
2. Main switch
3. Peristaltik pump 1 + 2

### 2.3.2 Installation and Alignment of DCU4-Tower

The DCU4 tower is provided with height adjustable feet for the alignment.

**Carefully level the DCU4 tower before you start it up.**

Figure 2-3: Feet of the DCU4 tower

1. Height adjustable feet (4 pieces)

Please proceed as follows in order to level the DCU4 tower:

- Loosen the counter nuts of the feet.
- Carefully level the DCU4 tower by turning the feet out or in.
- Fasten the counter nuts of the feet.

### 2.3.3 Signal and Power Supply of the DCU4 Tower

For the signal and power supply the DCU4 tower has got the connections described as follows.  
They are on the back side of the DCU4 tower.

Figure 2-4: Signal and power supply of the DCU4 tower (rear view)

1. Main supply
2. Power supply of the recirculation pump
3. Connections for digital signals and 24 volts power supply of the filtration unit
4. Connections for analog signals of the filtration unit
5. Ethernet connection
6. External 4–20 mA Signal

## 2.4 Filtration Unit

### 2.4.1 Structure of the Filtration System

The filtration unit mainly consists of the following single components:

Figure 2-5: Front view of the filtration unit

- |                         |                        |
|-------------------------|------------------------|
| 1. Retentate connection | 6. Recirculation bag   |
| 2. Recirculation pump   | 7. Air-filter          |
| 3. Pressure transducers | 8. Filter holder       |
| 4. Feed line            | 9. Permeate connection |
| 5. Feed connection      |                        |


Figure 2-6: Side view of the filtration unit

- |                                      |
|--------------------------------------|
| 10. Over pressure protection (Bag)   |
| 11. Temperature sensor               |
| 12. Load cell                        |
| 13. Flow meter in the retentate line |
| 14. Flow meter in the permeate line  |

## 2.4.2 Assembly and Alignment of Filtration Unit

The filtration unit is provided with height adjustable feet for the alignment.

**Carefully level the filtration unit before you start it.**

Figure 2-7: Feet of the filtration unit

### 1. Height adjustable feet (4 ... 6 pieces)

For the alignment of the filtration unit please proceed as follows:

- Loosen the counter nuts of the feet.
- Level the filtration unit by turning the feet in or out.
- Fasten the counter nuts of the feet.

## 2.4.3 Inspection of the filtration unit

Before you start the filtration system you have to carry out the following inspections:

- Inspection of the tubing of the complete system.
- Inspection of all connections and seals (if visible).
- Inspection of all electrical cable connections at the measuring sensors (flow meter).

#### 2.4.4 Signal and Power Supply of the Filtration Unit

The filtration unit is provided with respective connecting lines for the signal and power supply.  
(► „Connector pin assignment at the DCU4 tower“ see chapter 2.3.3)


Figure 2-8: Signal and power supply filtration unit

1. Main supply
2. Power supply of the recirculation pump
3. Connections for digital signals and 24 volts power supply of the filtration unit
4. Connections for analog signals of the filtration unit
5. Ethernet connection
6. External 4–20 mA Signal

#### 2.4.5 Mounting of the Recirculation Container (Bag as Recirculation Container)


The bag stand is assembled at a weight measuring system (e. g. load cell).  
The recirculation bag is hung up at the bag stand.

**Load cells are very precise weighing devices. For the mounting and the handling of the load cells please observe:**

- The load cell must not be exposed to torsion.
- The maximum bearing force of the load cell must not be exceeded during mounting and operation.
- Incorrect handling and/or mounting may lead to damage of the load cell.

The bag stand is screwed directly at the weighing device (load cell).

1. Weighing device (load cell)
2. Bag stand
3. Recirculation Bag

Figure 2-9: Bag stand for the recirculation bag

## Bag Stand installation

The weighing device is equipped with a welded fastener. The bag stand is directly screwed over it at the weighing device (load cell).

For the assembly of the bag stand please proceed as follows:

- Apply the bag stand with the fastener at the screw-down surface of the load cell.
- Screw down the bag stand with the provided socket screws (4 pieces) at the load cell.

## 2.4.6 Installation of the presterilized Single use Bag Loop

When unpacking the components please check whether the packaging is intact and no damage to the individual components is visible. Also, make sure that during the installation no parts are being damaged or contaminated.

### Procedure

1. Installation of Self Contained Unit following the instructions listed below.
2. Hang the recirculation Bag on the Bag holder and mount pressure sensor on Bag holder.
3. Install the Pump element in the peristaltic pump.
4. Mount the Feed pressure sensor.
5. Connect Feed line to filter element using OPTA.
6. Connect Retentate line to filter element using OPTA.
7. Mount the Retentate pressure sensor.
8. Mount the Retentate valve body.
9. Mount the Retentate flow sensor.
10. Connect Permeate line to filter element using OPTA.
11. Mount the Permeate pressure sensor.
12. Mount the Permeate valve body.
13. Mount the Permeate flow sensor.

Further instructions see below

### Pressure sensor

1. Remove the membrane Protection element from the locator on the pressure sensor.
2. Mount the open locator together with the pressure dome on the pressure transducer. It is important to ensure that the locator is aligned as shown in Figure 1. By applying light pressure on the two bars the gap is expanded and the pressure dome reaches into the groove of the transducer for fixation.
3. Turn the locator to 90 degrees – shown in Figure 2

Figure 1

Figure 2

### Electronic positioning valve – Valve body installation

1. During assembly of the valve body, the electronic positioning valve must be in the position "100%".
2. Mount the Valve body to the connector (see Figure 3).
3. Fasten the Valve with Tri-clamp (see Figure 4).
4. The connection in between the actuator and the valve body is established by moving the actuator into the position "0%".

Figure 3

Figure 4

### Flow sensor

1. The sensor mounted on the flow cell as shown in Figure 5.

Figure 5

### Installation of the Pump element

1. Lay the double y-tube element into the open pump head and place the branches over the fixing screws.
2. Close the lid. It might be helpful to use a wrench (SW10).

Figure 6

Before the assembly and cleaning of the filter holder you must carefully read the operating manual of the Sartocon Slice Self Contained Unit! In the case of incorrect handling of the filter cassettes severe injuries may result or the filter cassette may be damaged.

You must only use the provided torque wrench or the hydraulic assembly to clamp the Self Contained Unit!  
- Incorrect handling may result cassette damage.

Keep the stud bolts and screw nuts clean and free of grease!  
Greased or polluted threads lead to wrong clamp forces which may destroy the cassettes.

For the installation or change of the filter cassettes follow the instructions in „Sartocon Slice Self Contained Unit “<sup>[1]</sup>, see bibliography at the end of the chapter.

1. Install the Self Contained unit in the filter holder and clamp it with the recommended torque, or with the appropriate hydraulic pressure.


Figure 7

- Tighten the Cassette with the recommended clamping force according to the selected clamping device.
- When installing the Cassette make sure that the grooves are in correct orientation (to accommodate the upper and lower rod sites).
- Keep the threads of the rods clean (avoid grease).

Do not forget hydraulic clamping the upper locking bars | rods

Clamping force for Sartocon® Slice Self Contained  
Sartocon® Slice Self Contained Holder

	Clamping force Operating parameters	Clamping force Operating parameters Torque wrench
Sartocon® Slice Self Contained	14–17 kN	20 Nm


## OPTA-Connector


1. Connect the female and male OPTA connectors with the same colour tag. These colour codes ensure the correct allocation of the ports.

Orange: Feed  
Blue: Retentate  
Yellow: Permeate

### Making the Connection


Male connector


Female connector

#### Step 1

Remove the protective caps by pulling on the flaps that are located at the bottom of each cap. After removal of the caps membrane tapes will be released. These membrane tapes should always hang down to assure the correct assembly of the male and female connector bodies.



Male connector

Female connector

#### Step 2


Retract the collar that is fixed on the female connector to prepare for assembly. Assemble the two connectors by sliding the two connector heads together until the two pins of the male connector click into place. Once the female and the male connectors are assembled they should not be disassembled. The integrity of the male and female connectors is ensured by the membrane sterile barrier. Each membrane tape is equipped with a mating gripper to assure simultaneous removal of the sterile barrier.


### Step 3

After the two pins of the male connector are "clicked" into place, join the grippers at the bottom of each membrane tape together. The grippers should also "click" into place. Remove the membrane horizontally.

### Step 4

(A) Slide the collar that has been retracted back in place and (B) turn it clockwise until the lock "clicks" into place.

**Once the cassette and the assembly is connected we recommend to flush the loop.**

### Flushing

- Fill the recirculation bag with 5l of sterile water (WFI). Retentate and Permeate valves are fully open. Set delta P to a value of 1 bar. Flushing the system to the drain.

Refill the recirculation bag with 10l of sterile water (WFI) and pause for 15 min. Set delta P to a value of 1 bar.

- Restrict the retentate to an extent that a ratio between the Permeate flux to Retentate of 1:1:

Discard the permeate and retentate

3. Operation of the Filtration System

3.1 Starting the Unit | Emergency Stop Function

3. Operation of the Filtration System


Figure 3-1: Main switch

**Before you start the unit you have to read chapter 2 „Installation and Starting-up“ carefully.**

- For the starting of the unit you must proceed as follows:
1. Turn the main switch (pos. 1) to position "1" (vertical switch position).
  2. Await the system start (starting-up of the unit).

**In the case of eventual alarms you have to proceed as described in chapter 3.3.3 "Menu Bar → „Alarms“.**

3.2 Display

**For the operation of the touch panels you must not use pointy items (e. g. pens etc.).  
For the cleaning of the touch panel you must only use suited cleaning agents.**

3.2.1 Touch Fields of the Start Screen

After starting the unit the start screen (entry screen) shows the following windows:


Figure 3-2: Start screen of the SARTOFLOW® Alpha plus SU

The start screen is divided into the following frames:

1. Menu bar: The menu bar serves for the selection of the main menus (main, trend, settings, remote, alarm, calibration etc.).
2. Main window: Graphic of the process including all process values as well as the parameters of the control cycles.

### 3.3 Modes of Operation



#### 3.3.1 Manual Operation

This mode of operation allows a free configuration of all parameters (e. g. controller settings, alarms etc.) considering the permissible performance data (see chapter 4.1).

#### 3.3.2 Operation

It leads the user step by step through the configuration of the parameter settings needed for the selected sequence.

#### 3.3.3 Menu Bar


By pressing the button „Main“ the display in the main window changes to the display of the process. Here all process values as well as the control and alarm parameters are indicated.

The „Main“ menu can always be accessed from all menus and submenus.

The main display represents the active process. The indicated process values can be read off completely.

All settings relevant to the process (controller settings, process values) are controlled via the main display.


## Setting of the controller parameters

By pressing the button the display in the lower help window changes to the display where the controller values are entered.

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
Enter the alarm low and high limits

## Setting of the process values

By pressing the button the display in the help window at the bottom changes to the display where the process limit values are entered.

Enter the alarm high limit

Actual value display

Enter the alarm low limit

Switch on | off alarm


## Calibration of peristaltic pumps – P3000 & P3100


By pressing the button "P3000\_T", the display changes

Press the 'off' Button the display changes.


Insert tubing in the peristaltic pump. Use a glass cylinder to measure the pumped volume.

The calibration starts after you press 'CALIBTATE' in the next window.

By pressing OK the calibration stops.



Enter the measured the volume.


Confirm by pressing OK.

## Calibration of flow meter

By pressing the button "Calibration", the display changes and the phase for flow meter calibration starts.


Calibration values are displayed for Retentate and Permeate.

Enter the calibration factor k that is delivered with the disposable part.

The calibration factor k for the installed turbine is adjusted and confirmed by pressing OK.


## Tare function and Phases

By pressing the button Tare Pressure or Weight, the display changes and the phase for zero point tarring starts.

Further information Chapter: Phases

## Trend

By pressing the „Trend“ button the display in the main window changes to the display of the trend. Here all process values are represented as a chronological sequence.

The „Trend“ menu can always be accessed from all menus and submenus.

Settings at controller parameters and process values are not possible at this representation.


Figure 3-4: Trend display of the SARTOFLOW® Alpha plus SU


1. History (time)
2. Display of the process value size depending on the set scale
3. Select button (configurable)

## Configuration of select button

By pressing the button an additional window appears in the trend display. Here you find characteristics which can be assigned to the button after the respective selection.



Selection of the process value to be represented  
Selection of the colour of the trend line

Determine scale of the process value

Determine scale of the process values to the factory setting

## Settings

All system parameters, system settings etc. can be modified at the menu point "Settings". After pressing the button "Settings" the submenu points are displayed on the main window.

**Any settings and modifications must only be carried out by instructed and authorized persons as well as by the service personnel of Sartorius Stedim Systems GmbH.**

**Incorrect settings may cause errors or damages of the unit and may lead to eventual resulting accidents.**


Figure 3-5: Trend display of SARTOFLOW® Alpha plus SU

1. System parameters
2. PV ranges
3. External
4. Manual Operation
5. Trend
6. Service

Settings of the system parameters (e. g. time, date, IP address etc.)  
By pressing the button „system parameters“ the display in the main window changes to the display where the values are entered.

**System Parameters** [X]

Time  
12:50:10

Date  
2010-01-08

Failtime  
2:00:00

Screensaver  
0:00:00

Internet Config  
192.168.137.200

Settings of the. time, date, IP address etc.

Settings of the process values (e. g. min., max., alarm values etc.)  
By pressing the button „process values“ the display in the main window changes to the display where the values are entered.

**Process Value Ranges** [X]

Ch.	Process Value	Min	Max	Alarm	Alarm Low	Alarm High
1	TIR2100	0.0 °C	150.0 °C	enabled	0.0 °C	150.0 °C
2	FIRC2600	0.00 bar	5.00 bar	disabled	0.00 bar	3.00 bar
3	FIRC2500	0.00 bar	5.00 bar	disabled	0.00 bar	3.00 bar
4	FIRC2700	0.00 bar	5.00 bar	disabled	0.00 bar	3.00 bar
5	FIRC2600	0.0 L/h	1500.0 L/h	disabled	0.0 L/h	1500.0 L/h
6	FIRC2700	0.0 L/h	150.0 L/h	disabled	0.0 L/h	120.0 L/h
7	WIRC2100	-5.00 kg	15.00 kg	enabled	-5.00 kg	10.10 kg
8	P3000	0.0 %	100.0 %	disabled	0.0 %	100.0 %
9	P3100	0.0 %	100.0 %	disabled	0.0 %	100.0 %
10	TMP	0.0 bar	4.0 bar	disabled	0.0 bar	4.0 bar
11	DPRESS	0.0 bar	5.0 bar	disabled	0.0 bar	5.0 bar
12	QIR2702	0.000 CU	5.000 CU	disabled	0.000 CU	5.000 CU
13	P3000_T	0.00 L	500.00 L	disabled	0.00 L	500.00 L

Settings of the process values (e.g. min., max., alarm values etc.)

Settings of external connections (not available)

Configuration of all analog and digital entries and exits as well as control cycles. By pressing the button „manual operation“ the display in the main window changes to the display where the values are entered and all process values as well as in- and outputs are indicated.

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User settings.


Edit User.

User	1
Real Name	1
Change PW	
Group	level1
Dis/Enable	enabled
Expire	1970-01-01
Delete	ok

Log Book.


Log Book with process messages.  
Changes and operations are recorded with time and user.



Service menu  
This menu can only be operated by Sartorius Service.


## Alarm

By pressing the „Alarm“ button the display in the main window changes to the display of alarms. Here all alarms which have occurred during the running process are displayed.

Status „Ready for operation“

Status „Alarm“

Alarm

12345

Time	Message	State	Ack.	Reset
2010-01-08 12:51:55	Process Value WIRC2100 Low Alarm At 0.05 kg Alarmlimit: Min 1.00 Max 10.10	UNACK	ACK	RST
2010-01-08 12:51:37	Process Value TOTAL Low Alarm At 0.00 L Alarmlimit: Min 5.00 Max 1000.00	UNACK	ACK	RST



ACK ALL6

Figure 3-6: Alarm display of the SARTOFLOW® Alpha plus SU


1. Alarm time (time of the arising alarm).
2. Alarm text (description of the alarm).
3. Alarm status (unreceipted | receipted).
4. Button to receipt an alarm (single receipt).
5. Button to cancel an alarm (receipt and cancel).
6. Button for the actualization of the alarm display.

## Remote mode

The remote mode is deactivated. There is a communication between the external control and the DCU4 control unit. The remote control is deactivated. The operator cannot access the DCU4 control unit from an external control (e. g. MFCS/win). Data transfer to an external control will not be interrupted.

The remote mode is activated. There is a communication between the external control and the DCU4 control unit. The bi-directional control is activated. A remote control of the DCU4 control unit is possible.

If you use MFCS/win or MFCS/DA the remote mode should be activated, in order to set the DCU4 control unit in the bi-directional status. Only that way it is guaranteed that the remote control of the DCU4 control unit works correctly.

## Phases


Before starting a process you should tare the pressure sensors and the load cell. The flow meters have to be calibrated by entering the k factor in the calibration window.

## General

The system provides the following phases and controller:

- |                  |                  |           |
|------------------|------------------|-----------|
| 1. Tare Pressure | 4. Concentration | 7. DPRESS |
| 2. Tare Wight    | 5. Diafiltration |           |
| 3. Conditioning  | 6. TMP           |           |


The phases will be started from the operator panel and the DCU executes the steps with interaction of the operator. A phase can be started and stopped at each state by the operator. If no phase is active, the system is in manual mode. In this mode the feed pump can be started and all valves can be operated directly from the operator panel.

### 1. Phase „TARE PRESSURE“

The sequence “PRESSURE” is used to set all the pressure sensors to zero automatically before using another phase. It has to be started after a new bag assembly is installed. This phase should only be started with an empty bag immediately after the installation. The sequence executes the steps automatically, also triggers “operator information | messages or questions” and waits for operator interaction wherever it is necessary as described below: This phase starts after activating the start button in the phase menu and confirming by pressing OK.


Sequential Function Chart

### 2. Phase „TARE WEIGHT“

The sequence “WEIGHT” is used to set the load cell to zero automatically before using another phase. The sequence executes the steps automatically, also triggers “operator information | messages or questions” and waits for operator interaction wherever it is necessary as described below: This phase starts after activating the start button in the phase menu and confirming by pressing OK.

Sequential Function Chart


Parameter table

### 3. Phase „CONDITIONING“


The phase "Conditioning" is used to flush the cassettes with buffer and prepare them for the process. The peristaltic pumps fill the recirculation bag with buffer until the filling volume reaches the defined end level. The buffer flushes the filter cassettes and is collected in a waste container. The sequence executes the steps automatically, also triggers "operator information | messages or questions" and waits for operator interaction wherever it is necessary. This phase starts after activating the start button in the phase menu and confirming by pressing OK. The phase should only be used for the flushing new bag assemblies after the installation.

## 38


### Operation of the Filtration System

## 5. Phase „DIAFILTRATION“


The sequence "Diafiltration" substitutes media for buffer. During the concentration phase the working volume decreases to a defined level. Once the hold level is reached the peristaltic pump will adjust the weight to the defined level. If the level sensor registers a weight that is 0,1 kg below the defined level, the peristaltic pump stops and the DCU submits a query to the operator and asks for more buffer. When the buffer supply can be guaranteed, the operator has to confirm the status and the system fills up the tank until the set point is reached and continues with the diafiltration step. The sequence executes the steps automatically, also triggers "operator information | messages or questions" and waits for operator interaction wherever it is necessary as described below: This phase starts after activating the start button in the phase menu and confirming by pressing OK.


6. „TMP CONTROL“

The sequence “TMP Control” is used to control the transmembrane pressure by using the “difference pressure control loop” and the “TMP control loop”. The sequence executes the steps automatically, also triggers “operator information | messages or questions” and waits for operator interaction wherever it is necessary as described below: This phase starts after setting the TMP controller into “manual” or “automatic” mode.

7. „DPRESS CONTROL“

The sequence “DPRESS” is used to control the differential pressure by using the “DPRESS control loop”. The sequence executes the steps automatically, also triggers “operator information | messages or questions” and waits for operator interaction wherever it is necessary as described below: This phase starts after setting the DPRESS controller into “manual” or “automatic” mode.

For your notices:

---

---

---

---

---

## 4. Appendix

### 4.1 Technical Data

## 4. Appendix

Dimensions	1300 × 1200 × 960 mm
Weight	Approx. 150 kg
Power supply	110 Volt AC or 230 Volt AC (16 A)
Pressure	–1 to +3 bar
Filter area	0.1 m <sup>2</sup> to 0.3 m <sup>2</sup>
Weighing range	0 to 100 kg
Pump capacity	300 L/h at 3 bar
Material	Stainless steel
Communication connection	Ethernet

### 4.2 Cleaning and Care

The cleaning and maintenance intervals mainly depend on the degree of wear and pollution of the accessories.

You can clean the machine surfaces with common household dish liquids or alcohol.

For the care of stainless steel surfaces the respective stainless steel cleaning agents of leading manufacturers are suited.

**Please avoid aggressive cleaning agents, e.g. chloric cleaning agents**

### 4.3 Maintenance and Service

The maintenance of the SARTOFLOW® Alpha plus SU is restricted to the general test and inspection works of the filtration module.

► See chapter 2.4.3 „Inspection of the Filtration Unit“

In the case of other malfunctions of the unit please contact your responsible service of Sartorius Stedim Systems GmbH.

► See chapter 4.10 „Contact – Addresses“

#### 4.4 Causes of Malfunction

Malfunction	Possible reasons	Possible countermeasures
Display rests dark	Unit has not been switched on	Check power supply   start the unit again
	Incorrect system start	Start the system again
	Emergency stop is activated	Deactivate emergency stop
Unit cannot be started	Wrong power supply	Check power supply
	Emergency stop is activated	Deactivate emergency stop
Pump P-2500 does not start	Pump cover is not closed	Close the pump cover
Flexible-tube pump P-3000 does not start	Incorrect controller settings	Check controller settings
Sequence cannot be started	Other sequence is active	Check active sequences and eventually deactivate them
	Alarms are active	Receipt alarms
Measured values are wrong or they are not displayed	Connection error between measuring device and DCU4 control unit	Check cable connections
No communication between DCU4 control unit and MFCS	Remote button has not been activated	Check remote settings
	Ethernet cable has not been connected correctly	Connect Ethernet cable correctly
Measured values indicate the maximum measured value	Cable break	Check cable connections
	Measurement range has been exceeded	Check measurement range
Wrong results at pressure measurement	Measuring instrument not calibrated	Calibrate measuring device
	No zero point alignment	Carry out zero point alignment
Wrong results at level measuring	Weighing machine not tared	Tare weighing machine
	Measuring instrument not calibrated	Calibrate measuring instrument
Wrong results at flow measuring	The calibration factor k has not been given correctly	Adjust the calibration factor k to the installed turbine

### Disassembly

► „Fill out the declaration of decontamination (chapter 4.8)“

Pay attention to the qualification of the disassembling personnel.

- Only instructed | skilled personnel (experts) must demount the filtration system.
- Only qualified personnel must be in charge with specific areas (hydraulic systems, gas, electro technology etc.).

Conditions for the disassembly.

- The filtration system must be completely emptied before disassembling it.
- Strong pollutions must be removed. You must clean the components polluted by dangerous substances (before disassembling).
- Set the main switch to „0“ and save this position.
- Lock the energy supply and save it.
- Separate the DCU4 tower from the power supply unit.

Unload the containers and lines (under pressure).

- Unload the containers and lines under pressure at the specified locations.
- Close the connections after having disconnected the lines.

Springs and similar components (under tension) must be unloaded.

- Unload the tightened springs.

Carry out the disassembly.

- Arrange the disassembled filtration unit components so that they are suited for the transport.

### Disposal

Collect dangerous substances.

- Operating materials may contain dangerous substances. Clean the components which have been polluted by dangerous substances (before the disassembly).
- Collect all operating materials and other harmful substances in collecting tanks and dispose them according to the statutory provisions (professional disposal).

Please observe the statutory provisions.

- The disposal of the filtration system must be carried out according to the applicable laws and regulations.
- Since the filtration system is fixed in the operating state, the old electronics guideline „WEEE“ is not applicable.
- You must eventually sign off the filtration unit or any parts must | can be returned to the manufacturer.
- Before the dispatch you must identify the parts which have been polluted by dangerous substances.

You have to care for a separate disposal of the different components according to the type of material.

- Metal and non-ferrous metals at metal recycling companies.
- Plastics and compounds at plastic recycling companies.

#### 4.6 Warranty

As far as not stipulated otherwise all products of Sartorius Stedim Systems GmbH are subject to a warranty according to the valid general terms of business.

- All components of SARTOFLOW® Alpha plus SU have been designed for the use in common filtration processes. Under special environmental ambient conditions or when using aggressive media you must check the suitability and durability of all components.
- The warranty applies for all construction, fabrication or material faults and the resulting malfunctions. It includes the repair or the replacement of damaged parts.
- The warranty does not refer to consumables and parts which are subject to wear and tear (for example electrodes, o-rings, seals, membrane filters). Excluded from the guarantee are also any malfunctions, defects or resulting damages due to corrosion at improper ambient conditions or at the use of corroding materials.

Warranty and liability claims are excluded at those components,

- which do not correspond with the specifications for the unit and its application.
- which are procured from other suppliers.
- for which Sartorius Stedim Systems has not given its written authorization of use.

#### 4.7 EC Declaration of Conformity

EC declaration of conformity see enclosure

#### 4.8 Declaration of Decontamination

Declaration of decontamination see enclosure

#### 4.9 Contact Addresses

Contact addresses see enclosure

#### 4.10 Information and Instructions on Disposal

Attention must be paid to the applicable laws, rules and regulations, when disposing of the contaminated components.

The packaging is to be taken to a local waste disposal site if no longer required. The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste. EU legislation requires its Member States to collect electrical and electronic equipment and disposed of it separately from other unsorted municipal waste with the aim of recycling it.

In Germany and several other countries, Sartorius Stedim Biotech GmbH itself assumes responsibility for the return and conformant disposal of its electronic and electrical products. These products may not be placed with household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other member nations of the European Economic Area (EEA), please contact our local service technicians or our Service Center in Goettingen, Germany:

Sartorius Stedim Biotech GmbH  
Servicezentrum Biotechnologie  
August-Spindler-Strasse 11  
37079 Goettingen, Germany

WEEE-Reg.-No. DE 89907997

In countries that are not members of the European Economic Area (EEA) or where no Sartorius Stedim Biotech subsidiaries or dealerships are located, please contact your local authorities or a commercial disposal operator.

Prior to disposal and |or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.

Sartorius Stedim Biotech will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the package leaflet | the included CD ROM or to our web site ([www.sartorius-stedim.com](http://www.sartorius-stedim.com)) for more detailed information and addresses for repair service or disposal of your device.

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NEW USER    GROUPS    PARAMETERS

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
Edit User 1	
User	1
Real Name	1
Change PW	
Group	level1
Dis/Enable	enabled
Expire	1970-01-01
Delete	ok

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



A screenshot of the bottom navigation bar of the application. It contains seven icons in a row: 'Main' (a house icon), 'Trend' (a line graph icon), 'Settings' (a gear icon), 'Remote' (a remote control icon), 'Alarm' (a bell icon), 'User' (a person icon), and 'Shutdown' (a power button icon).

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
## Phasen

## Allgemeines

1. Tare Pressure	4. Concentration	7. DPRESS
2. Tare Weight	5. Diafiltration	
3. Conditioning	6. TMP	

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Copyright by  
Sartorius Stedim Biotech GmbH,  
Goettingen, Germany.  
All rights reserved. No part of  
this publication may be reprinted  
or translated in any form or by  
any means without the prior  
written permission of Sartorius  
Stedim Biotech GmbH.  
The status of the information,  
specifications and illustrations  
in this manual is indicated by  
the date given below.  
Sartorius Stedim Biotech GmbH  
reserves the right to make  
changes to the technology, fea-  
tures, specifications and design  
of the equipment without notice.

Status:  
October 2012,  
Sartorius Stedim Biotech GmbH,  
Goettingen, Germany